
ORTHOKIT 1.0

EN · DE · NL · FR · ES


ORTHOKIT 1.0

CONTENTS

ORTHOKIT 1.0 EN	3
ORTHOKIT 1.0 DE	12
ORTHOKIT 1.0 NL	22
ORTHOKIT 1.0 FR	32
ORTHOKIT 1.0 ES	42

ORTHOKIT 1.0 EN

• Eggbar Shoe	3
• Heartbar Shoe	4
• Straightbar Shoe	5
• Straightbar Performance Shoe	6
• Wide Toe Shoe	7
• Wide Branch Shoe	8
Dimensions & Sizes	9
Comments / Disclaimer	10
	11

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
WERKMAN ORTHOKIT 1.0

June 2013 saw the launch of a comprehensive scientific study on the subject of shoeing, funded by the Werkman Hoofcare Company and led by Dr. Jenny Hagen (Veterinary Anatomical Institute (Veterinär-Anatomisches Institut) of Leipzig University).

The study focused on the impact of various horseshoes and shoeing modifications, taking into account diverse soil types, on:

- the distribution pattern of the vertical compressive forces exerted on the hoof;
- the size of the load-bearing area;
- the gait pattern;
- the orientation of the phalanges.


The results of this study were used for the development of the Werkman ORTHOKIT.

The Werkman ORTHOKIT provides you with an effective and straightforward means to turn Werkman horseshoes into various therapeutic shoes.

This brochure explains how the ORTHOKIT modifications are made and how they function.

EGGBAR SHOE

When using the bar on Warrior shoes with side or quarter clips or on the Warrior Special shoe, the effect of the modification is enhanced by the shape of the toe of the shoe.

FIELD OF APPLICATION

- Pathologies of the deep flexor tendon
- Pathologies of the check ligament of the deep flexor tendon
- Chronic navicular disease
- Palmar foot pain
- Inflammations of the ligaments of the toe


EFFECT

- Expansion of the palmar support area and of the palmar lever arm
- Increase of the angle of the coffin bone on soft soil
- Reduction of the leverage exerted on the front section of the (toe of the) hoof on all soils
- Reduction of the tension of the deep flexor tendon and a reduced load on the navicular area
- Without frog contact: protection of frog and palmar hoof section (navicular bone, attachment of the deep flexor tendon) on all types of soil
- With frog contact: continual increase of the load-bearing surface in the palmar hoof section, regardless of the type of soil

This Eggbar can easily be modified into an open toe shoe.

HEART BAR SHOE

When using the bar on Warrior shoes with side or quarter clips or on the Warrior Special shoe, the effect of the modification is enhanced by the shape of the toe of the shoe.

FIELD OF APPLICATION

- Founder/Laminitis
- Pathologies of the deep flexor tendon
- Pathologies of the check ligament of the deep flexor tendon
- Chronic navicular disease
- Palmar foot pain
- Inflammations of the ligaments of the toe
- Hoof cracks


EFFECT

- Expansion of the palmar support area without modification of the palmar lever arm
- Reduction of the risk of losing the horseshoe due to pulling the shoe (in comparison with shoeing with palmar extension)
- Increase of the angle of the coffin bone on soft soil
- Reduction of the tension of the deep flexor tendon and a reduced load on the navicular area
- Without frog contact: protection of frog and palmar hoof section (navicular bone, attachment of the deep flexor tendon) on all types of soil
- With frog contact: continual increase of the load-bearing surface in the palmar hoof section, regardless of the type of soil

This Heart Bar can easily be modified into an open toe shoe.

STRAIGHT BAR SHOE

When using the bar on Warrior shoes with side or quarter clips or on the Warrior Special shoe, the effect of the modification is enhanced by the shape of the toe of the shoe.

FIELD OF APPLICATION

- Founder/Laminitis
- Pathologies of the deep flexor tendon
- Pathologies of the check ligament of the deep flexor tendon
- Chronic navicular disease
- Palmar foot pain
- Inflammations of the ligaments of the toe
- Hoof cracks


EFFECT

- Expansion of the palmar support area without modification of the palmar lever arm
- Reduction of the risk of losing the horseshoe due to pulling the shoe (in comparison with shoeing with palmar extension)
- Increase of the angle of the coffin bone on soft soil
- Reduction of the tension of the deep flexor tendon and a reduced load on the navicular area
- Without frog contact: protection of frog and palmar hoof section (navicular bone, attachment of the deep flexor tendon) on all types of soil
- With frog contact: continual increase of the load-bearing surface in the palmar hoof section, regardless of the type of soil

This Straight Bar can easily be modified into an open toe shoe.

STRAIGHT BAR PERFORMANCE SHOE

When using the bar on Warrior shoes with side or quarter clips or on the Warrior Special shoe, the effect of the modification is enhanced by the shape of the toe of the shoe.

FIELD OF APPLICATION

- Founder/Laminitis
- Pathologies of the deep flexor tendon
- Pathologies of the check ligament of the deep flexor tendon
- Chronic navicular disease
- Palmar foot pain
- Inflammations of the ligaments of the toe
- Hoof cracks


EFFECT

- Expansion of the palmar support area without modification of the palmar lever arm
- The curved bar reduces the risk of losing the horseshoe due to pulling the shoe (in comparison with shoeing with palmar extension)
- Increase of the angle of the coffin bone on soft soil
- Reduction of the tension of the deep flexor tendon and a reduced load on the navicular area
- Without frog contact: protection of frog and palmar hoof section (navicular bone, attachment of the deep flexor tendon) on all types of soil
- With frog contact: continual increase of the load-bearing surface in the palmar hoof section, regardless of the type of soil

This Straight Bar Performance shoe can easily be modified into an open toe shoe.

WIDE TOE SHOE

By making the branches of the shoe narrower (grinding, fuller cutting) the effect of the horseshoe is enhanced on soft soil.


FIELD OF APPLICATION

- Pathologies of the superficial flexor tendon
- Pathologies of the check ligament of the superficial flexor tendon
- Pathologies of the suspensory ligament and its check ligaments

EFFECT

- Increase of ground contact area in the toe
- Reduction of the sinking in of the toe on soft soil
- Less upright orientation of the dorso-palmar coffin bone angle than with standard shoes on soft soil
- Position of the coffin bone matches the orientation on hard soil = constant
- On soft soil reduction of the tension of the superficial flexor tendon, the check ligament of the superficial flexor tendon and the suspensory ligament

WIDE BRANCH SHOE

The inner or outer horseshoe branch can be widened by welding in a wide branch bar.

By narrowing the opposite branch (grinding, fuller cutting), the effect of the shoe is enhanced on soft soil.

Grinding/forging the outside (ground side) of the narrower branch ("rolling") intensifies the effect on any type of surface. Both can be combined.

FIELD OF APPLICATION

- Damage to the collateral ligaments (widening on the side of the pathology / narrowing on the opposite side / rolling on the opposite side)
- Unilateral osteoarthritis of the toe joints (widening on the side opposite of the pathology / narrowing on the side of the pathology / rolling on the side of the pathology)


EFFECT

- One-sided increase of ground contact area
- Increased sinking in of the opposite side on soft soil
- Relief of collateral ligaments on the widened side
- Relief for unilateral osteoarthritis on the side that was not widened

DIMENSIONS

ORTHOKIT 1.0


Model	Size	mm	
		Height	Width
Ortho S	Wide Toe Shoe	68.5	17.2
	Wide Branch Shoe	89.5	21.4
	Heart Bar Shoe	62.2	75.9
	Straight Bar Shoe	61	33.9
	Straight Bar Performance	61	33.9
	Eggbar Shoe	70.7	36.8
Ortho M	Wide Toe Shoe	79.3	20.1
	Wide Branch Shoe	102.8	25.1
	Heart Bar Shoe	68.7	84.1
	Straight Bar Shoe	66.6	37
	Straight Bar Performance	68	38
	Eggbar Shoe	79.2	40.8
Ortho L	Wide Toe Shoe	90.2	22.5
	Wide Branch Shoe	118	28.4
	Heart Bar Shoe	78.6	96.9
	Straight Bar Shoe	75.4	42.6
	Straight Bar Performance	75.4	42.6
	Eggbar Shoe	89.7	46.5

EGG BAR SHOE


HEART BAR SHOE


STRAIGHT BAR SHOE

STRAIGHT BAR PERF. SHOE

WIDE TOE SHOE

WIDE BRANCH SHOE

COMMENTS

SIZES

Size S fits flush for a Warrior front shoe size 0
(Equates to Ryder Original Size 2 / Euroskill size 0)

Size M fits flush for a Warrior front shoe size 2
(Equates to Ryder Original Size 4 / Euroskill size 2)

Size L fits flush for a Warrior front shoe size 5
(Equates to Ryder Original Size 7 / Euroskill size 5)

IMPORTANT

The ORTHOKIT modifications have different effects on different types of soil. It is important to consider the bedding in the box and type of soils that the horse is worked on for optimum effect.

Keep in mind that relieving the load on an anatomical structure will always be accompanied by an increase in the load in another area.

The ends of the egg bar, heart bar, straight bar and performance bar are premade as left and rights. Make sure that the bars are dropped in the right way.

DISCLAIMER

The Werkman ORTHOKIT should only be selected, used, and mounted by farriers and veterinarians who possess adequate knowledge of orthopedic shoeing. Inexpert or erroneous selection or mounting of the Werkman ORTHOKIT (especially when welding the parts together) can undo the desired effect and even cause damage. Therefore, C.V. Werkman disclaims any liability for potential injuries, damage, costs or other consequences, caused by inexpert or erroneous selection and/or mounting of the Werkman ORTHOKIT.

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© WERKMAN HOOFCARE BV