





# ORTHOKIT 1.0 & 2.0

## ORTHOKIT 1.0


In June 2013, a comprehensive scientific study on the subject of shoeing was launched, funded by the Werkman Hoofcare company and led by Dr. Jenny Hagen (Veterinary-Anatomical Institute of Leipzig University). The results of this study were used for the development of the Werkman ORTHOKIT.

The Werkman ORTHOKIT provides you with an effective and straightforward means to turn Werkman horseshoes into various therapeutic shoes.


ORTHOKIT 1.0 is also available in ALUMINIUM.




Eggbar Shoe




Heart Bar Shoe




Straight Bar Shoe




This Eggbar, Heart Bar, Straight Bar and Performancebar can easily be modified into an open toe shoe.



Straight Bar Performance





Wide Toe Shoe




Wide Branch Shoe

# ORTHOKIT 2.0


The "ORTHOKIT 2.0" remains true to the original ideas, enlarging the range of easy to use orthopaedic horseshoe inserts designed by the experience of those in the field and supported by scientific research from the University of Leipzig.




Z-Bar Shoe




Central Bar Shoe



Spider Plate Shoe



Lateral Extension Shoe



Heart Bar Shoe



**WERKMAN HOOF CARE**

Helpermolenstraat 43  
NL 9721 BT Groningen

Postbus 6002  
NL 9702 HA Groningen

T +31 (0)50 52 18 121  
F +31 (0)50 52 58 256

shop@werkmanhoofcare.com  
werkmanhoofcare.com